
PROMESH[®] ABSOLIGHT **PROMESH ABSOLIGHT**

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Partially Absorbable Mesh – Poliglecaprone PGA-CL / Polypropylene

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Inguinal hernias by laparoscopy or open surgery and incisional hernias by open surgery

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Partially absorbable mesh – Innovative knitted fabric made of partially absorbable thread which offers greater strength with a lightweight implanted material

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PGA-CL (Poliglecaprone) and Polypropylene thread structure

- Bicomponent monofilament composed of alternation of non absorbable polypropylene and absorbable PGA-CL sections

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
Controlled absorption


- Surfacic mass before absorption: 60 g/m² (mesh thickness: 0.55 mm)
- Absorption time: 90-120 days
- Surfacic mesh post absorption: 28 g/m²
Minimal residual material is implanted in the body (mesh thickness: 0.4 mm)
- After absorption of the PGA-CL, the monofilament structure becomes a multifilament structure

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Optimal physical characteristics after absorption

- Mechanical resistance adapted to the pressures on the abdominal wall, in any stress axis^[1]
- Optimised tissue colonisation thanks to adequate pores size and low shrinkage rate of mesh after implantation, approx. 1.9%

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Easy implantation

- The excellent memory shape facilitates the deployment and the handling in extra peritoneal area, by laparoscopy or open surgery
- The colored stripes provide in a better visibility intraoperative procedure
- Sizes and geometries are adapted for Lichtenstein technique or laparoscopic procedures
- Compatible with usual fixation devices: glue, sutures, staples, helicoïdal fixations, absorbable or not

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Reduction of pain

- PGA-CL generates a moderate tissue reaction during the absorption
- The reduced postoperative pain facilitates faster rehabilitation of the patient and improves his quality of life^[3]

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Security and strength after absorption


- The knitted large pore fabric improves tissue integration of the implant thus forming a solid fibrous connective tissue^[2]
- The remaining structure of polypropylene multifilament provides a long-term effective reinforcement
- The flexible and lightweight prosthesis is compatible with the mobility of the abdominal wall and improves patient comfort^[4]

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Partially absorbable mesh
PGA-CL / Polypropylene

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Literatura: / Bibliography:

- [1] Etude RDE 13-54 (2014), Service R&D / RDE study 13-54 (2014), R&D department
 [2] Amid PK. Classification of biomaterials and their related complications in abdominal wall hernia surgery. *Hernia* 1997;1:951-960
 [3] Rapport HAS : avis de la Commission du 18 novembre 2008 / HAS report: opinion of the Committee dated 2008, November 18th
 [4] Note NRD 15-10 (2015), Service R&D / report 15-10 (2015), R&D department

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