

PROMESH[®] LIGHT **PROMESH LIGHT**

Dwuwymiarowa siatka lekka – Luźno pleciony polipropylen

Two-dimensional Mesh – Light Knitted Polypropylene

SIATKI PRZEPUKLINOWE / MESHES


Do zaopatrywania przepuklin metodą laparoskopową lub otwartą
Hernias by laparoscopy or open surgery

PROMESH® LIGHT



Dwuwymiarowa siatka lekka zbudowana z przeplatanych, monofilamentowych włókien polipropylenu:

- Gramatura: 30 g/m², wielkość porów: 1,5 x 1,5 mm
- Grubość siatki: 0,32 mm, grubość włókna: 0,08 mm
- Minimalna ilość implantowanego materiału w celu ograniczenia stanu zapalnego oraz redukcji bólu pooperacyjnego^[1]
- Wykonana z polipropylenu o udowodnionej tolerancji klinicznej^[2]



Light two-dimensional knitted fabric made of polypropylene monofilament:

- Surfacic mass: 30 g/m², pore size: 1.5 x 1.5 mm
- Mesh thickness: 0.32 mm, fiber thickness: 0.08 mm
- Low implanted material quantity to reduce the inflammatory reaction and the post-operative pain^[1]
- Polypropylene structure with proven clinical tolerance^[2]

Kod produktu Reference	Rozmiar [cm] Size [cm]	Kształt Shape
PPL0611	6,2 x 11,1	
PPL0913	9 x 13	
PPL1015	10,1 x 15,1	
PPL1515	15 x 15	
PPL1530	15 x 30	
PPL2030	20 x 30	
PPL3030	30 x 30	

Siatki PROMESH® Light nie mogą być docinane.
Skontaktuj się z nami w celu doboru odpowiedniego rozmiaru.
PROMESH® Light meshes must not be cut.
Please contact us for other sizes.

Literatura: / Bibliography:

[1] Brij B. Agarwal et al. Traditional polypropylene and lightweight meshes in totally extraperitoneal inguinal herniorrhaphy. International Journal of Surgery 8 (2010) 44-47.

[2] Yilmaz Bilisel and Ilker Abci. The search for ideal hernia repair; mesh materials and types. International Journal of Surgery (2012) 1-5.



MADENS MACIEJ MARCINKOWSKI, ul.Żeromskiego 9, 60-544 Poznań
www.madens.pl ; e-mail: biuro@madens.pl
Tel/fax: +48 61 843 48 00



Classe IIb - CE 0459 - LNE/G-MED - PETERS SURGICAL -

Ref.: PL01-98R31 / 07-2013