

EURODERM

EURODERM MC POWDER-FREE

HIGHLY SKIN-FRIENDLY -
the Semperit Post Production Process

EASY TO DON - patented synthetic inner coating

MAXIMUM COMFORT - fatigue-proof working



EURODERM MC POWDER-FREE



The EURODERM MC offers good reasons to switch to powder-free gloves.

HIGHLY SKIN-FRIENDLY

This powder-free glove provides improved safety for staff, patients and visitors as there are no powder-related dermatological, allergic and toxic health risks.

EASY TO DON

With its synthetic inner coating the glove is easy to put on due to its patented net structure.

SURE GRIP

The Euroderm MC has more grip and thus combines safety, tactical properties and instrument guidance.

FEELS LIKE A SECOND SKIN

Due to its resistance, the wall thickness of the Euroderm PF has been optimized to such a level that your tactile sensitivity is not in the least impaired, while at the same time providing maximum safety. The glove feels like a second skin and when double donning your sensitivity is not compromised by the addition of the second pair.

AVOIDING HAND FATIGUE

This relaxed hand posture should guarantee that you can work with the Euroderm MC without experiencing hand fatigue and still ensuring you are able to guide your instruments relaxed, calmly and safely, even in the longest of operations. The extrawide back of the hand also helps to achieve this.

TECHNICAL DATA

POWDER-FREE

TYPE	Sterile surgical glove for single use, powder-free, with synthetic lining			
MATERIAL	Natural rubber latex (NR)			
COLOUR	White			
GLOVE SHAPE	Fully anatomical with rolled rim			
OVERALL LENGTH/SIZE as per EN 455-2	5 ½, 6 and 6 ½	min. 265 mm		
	7, 7 ½ and 8	min. 275 mm		
	8 ½ and 9	min. 280 mm		
WALL THICKNESS measured in single layer in the palm area	0.19 – 0.23 mm			
IMPERMEABILITY / FORCE AT BREAK as per EN 455-1 / as per EN 455-2	AQL 0.65 / ≥ 12 N			
DURABILITY in original packaging if stored as per DIN 7716, ISO 2230	3 years			
STERILISATION	Radiation STERILE R			
PACKAGING	Left and right glove with turned up cuff in inner pouch, ozone-tight, sealed in flat medical peel pack.			
	In dispenser carton with sterilisation indicator: 50 pairs In transport carton with sterilisation indicator: 300 pairs			
ARTICLE NUMBERS	Size 5.5	822756524	Size 7.5	822756724
	Size 6.0	822756604	Size 8.0	822756804
	Size 6.5	822756624	Size 8.5	822756824
	Size 7.0	822756704	Size 9.0	822756904
LABELLING	EN 1041, EN 455, EN 420, EN 374, ISO 15223			

Semperit Investments Asia Pte Ltd,
Business Segment Sempermed
Singapore 609434 · 8 Jurong Town Hall Road
#29-03 to 06 The JTC Summit

Semperit Technische Produkte Gesellschaft m. b. H. Segment Sempermed
Austria · 1030 Vienna · Modecenterstraße 22
Tel.: +43 1 79 777-0 · Fax: +43 1 79 777-630
e-mail: sempermed@semperit.at · www.sempermed.com