


PROMESH® SOFT EVOLUTION PROMESH SOFT EVOLUTION

Trójwymiarowa siatka – Pleciony poliester
Three-dimensional Mesh – Knitted Polyester



Do zaopatrywania przepuklin metodą otwartą
Hernias and incisional hernias by open surgery


PROMESH® SOFT EVOLUTION

Dzianina 3D / 3D Fabric



Trójwymiarowa dzianina zbudowana z multiflamentowego poliestru o strukturze przestrzennych oczek w kształcie sześcianu:

- Gramatura: 125 g/m², wielkość porów: 1,2 x 1,4 mm
- Grubość siatki: 2,3 mm, grubość włókna: 0,07 mm
- Wysoka wytrzymałość mechaniczna
- Makroporowata, elastyczna oraz wysoce poręczna
- Idealnie dopasowana do struktur anatomicznych, by zwiększyć komfort codziennego funkcjonowania pacjentów



Three-dimensional knitted fabric made of hexagonal stitches of polyester multifilaments:

- Surfacic mass: 125 g/m², pore size: 1.2 x 1.4 mm
- Mesh thickness: 2.3 mm, fiber thickness: 0.07 mm
- High mechanical resistance
- Light and large pore structure, very flexible with « soft » handling
- Anatomical conformation and post operative comfort for the patient

Kod produktu Reference	Rozmiar [cm] Size [cm]	Kształt Shape
PEV0611FT	06 x 11	
PEV0812PC	08 x 12	
PEV0613NF	06 x 13	
PEV0913	09 x 13	
PEV0913PC	09 x 13	
PEV1015	10 x 15	
PEV1515	15 x 15	
PEV1520	15 x 20	
PEV3030	30 x 30	

Skontaktuj się z nami w celu dobrania odpowiedniego rozmiaru.
Please contact us for other sizes.

Literatura: / Bibliography:

[1] Bilse Y and Abci I. The search for ideal hernia repair; mesh materials and types. International Journal of Surgery 10(2012)317-321



Madens Maciej Marcinkowski, ul. Żeromskiego 9, 60-544 Poznań
www.madens.pl ; e-mail: biuro@madens.pl
tel./fax: 61 843 48 00



Classe IIb - 0459 - LNE/G-MED - PETERS SURGICAL -

Tel.: (+48) 22 462-42-52 – Fax: (+48) 22 867-57-70

Ref.: PL01-98R92 / 07-2014